Automatic Upper Arm Blood Pressure Monitor M7 Intelli IT (ALRU) Instruction Manual

Table of Contents

Introduction	EN2
Operational principle of work	EN2
Safety Instructions	EN2
Intended Use	EN2
Application	EN3
Intended Users	EN3
Circumstances When to Consult a Physician	EN3
Indication for Use	EN4
Contraindications for Use	EN4
Side Effects	EN4
Minor Repairs	EN4
Before the First Use	EN4
Important Safety Information	EN5
1. Know Your Monitor	EN9
1.1 Contents/Product Includes	EN9
1.2 Monitor and Arm Cuff	EN9
1.3 Display and Symbols	EN12
2. Preparation	EN18
2.1 Installing Batteries	EN18
2.2 Using the AC Adapter	EN19
2.3 Setting the Date and Time	EN20
2.4 Blood Pressure Measurement Tips	EN22
2.5 Applying the Arm Cuff	EN23
2.6 Sitting Correctly	EN25
2.7 Selecting a User ID	EN25
3. Using Your Monitor without a Smart Device	EN26
3.1 Taking a Measurement	EN26
4. Using Your Monitor with a Smart Device	EN31
4.1 Pairing Your Monitor with a Smart Device	EN31
4.2 Taking a Measurement	EN32
4.3 Disabling/Enabling Bluetooth	EN33
5. Using Memory Function	EN34
5.1 Viewing the Readings Stored in Memory	EN34
5.2 Viewing the Average Reading	EN35
5.3 Viewing Weekly Averages	EN36
6. Other Settings	EN38
6.1 Deleting All Readings Stored in Memory	EN38
6.2 Restoring Your Monitor to the Default Settings	EN39
7. Error Messages and Troubleshooting	EN41
8. Maintenance	EN48
8.1 Maintenance	EN48
8.2 Storage	EN48
8.3 Cleaning	EN49
8.4 Calibration and Service	EN49
9. Description of accessories, medical devices, or non-medical devices intended for use in combination with the automatic blood pressure monitor OMRON M7 Intelli IT.	EN50
10. Specifications	EN51
11. Guidance and Manufacturer's Declaration	EN55
12. Useful Information	EN66

Introduction

Thank you for purchasing the OMRON M7 Intelli IT (ALRU) Automatic Upper Arm Blood Pressure Monitor.

Operational principle of work

Oscillometric method:

The beat in the pulsation generated by the contraction of the heart is captured as the pressure inside the cuff to measure the blood pressure. If the cuff wrapped around the upper arm is pressurized sufficiently, the blood flow stops, but the beat of the pulsation is present and the pressure inside the cuff receives this and oscillates. Next, as the pressure inside the cuff gradually decreases, the oscillation of the pressure within the cuff gradually increases and reaches a peak. As the pressure within the cuff decreases further, the oscillation decreases from its peak. The pressure within the cuff and the relationship with the increase and decrease of the oscillation within the cuff in this series of processes are stored into memory, calculations are carried out, and the blood pressure value is determined. The oscillometric method does not determine the blood pressure value instantly like a microphone type automatic blood pressure gauge with the auscultation method, but rather determines it from the series of change curves as explained above. Therefore, it is not easily affected by external noise, an electric scalpel or other electro surgical instruments.

EN

Safety Instructions

This instruction manual provides you with important information about the OMRON M7 Intelli IT (ALRU) Automatic Upper Arm Blood Pressure Monitor. To ensure the safe and proper use of this monitor, READ and UNDERSTAND all of these instructions. **If you do not understand these instructions or have any questions, contact your OMRON retail outlet or distributor before attempting to use this monitor. For specific information about your own blood pressure, consult with your physician.**

Intended Use

This device is a digital monitor intended for use in measuring blood pressure and pulse rate in adult patient population. This monitor detects irregular heartbeats during measurement and displays corresponding indication along with measurement results.

EN2

The device can detect an irregular pulse suggestive of Atrial Fibrillation (Afib). Please note that the device is not intended to diagnose Afib. A diagnosis of Afib can only be confirmed by an Electrocardiogram (ECG). If the Afib symbol appears, consult with your physician.

Application

It is mainly designed for general household use.

Intended Users

This device is designed to measure blood pressure in adults 18 years and older, in which the circumference of the shoulder corresponds to the size of the cuff.

Circumstances When to Consult a Physician

OMRON recommends to always consult your physician before using home healthcare equipment.

Self-diagnosis of measurement results and self-treatment without professional medical supervision can lead to incomplete or incorrect diagnosis and/or under or over medication. In cases where the measurement results or treatment effect does not agree with expected results, OMRON always recommends to consult your physician.

Consult with your physician before using this monitor if you have common arrhythmias such as atrial or ventricular premature beats or atrial fibrillation; arterial sclerosis; poor perfusion; diabetes; pregnancy; pre-eclampsia or renal disease.

Consult with your physician before using this monitor on an arm where intravascular access or therapy, or an arterio-venous (A-V) shunt, is present because of temporary interference to blood flow which could result in injury.

Consult with your physician before using this monitor if you have severe blood flow problems or blood disorders as cuff inflation can cause bruising.

Consult your physician before using this monitor if you had a mastectomy.

If the irregular heartbeat continues to appear, we recommend you to consult with and follow the directions of your physician.

Indication for Use

This product can be used for monitoring blood pressure and pulse rate in order to confirm e.g. medication efficacy or lifestyle improvement programs. The product can be used for trending of these parameters in general home situation.

Contraindications for Use

The product should not be used in case the cuff location is damaged by trauma or if other treatment interferes with proper use. This product should not be used on patients on which it causes high anxiety as the results will not reflect the actual user's blood pressure values.

DO NOT apply the arm cuff on your arm while on an intravenous drip or blood transfusion.

DO NOT use this monitor on infants, toddlers, children or persons who cannot express themselves.

Side Effects

- DO NOT take measurements more often than necessary because bruising, due to blood flow interference, may occur.
- Inflating to a higher pressure than necessary may result in bruising of the arm where the cuff is applied.
- Stop using this monitor and consult with your physician if you experience skin irritation or discomfort.

Minor Repairs

Besides replacement of depleted batteries (where applicable) and re-fitting with approved accessories (mentioned in Instruction manual) due to wear and tear, no further minor repairs are to be carried out on the main product. If a defect occurs, consult your OMRON retail outlet or distributor as mentioned on the packaging.

Before the First Use

Remove this monitor from the packaging and inspect for damage. If this monitor is damaged, DO NOT USE and consult with your OMRON retail outlet or distributor.

Important Safety Information

Read the Important Safety Information in this instruction manual before using this monitor.

Follow this instruction manual thoroughly for your safety.

Keep for future reference. For specific information about your own blood pressure, **CONSULT WITH YOUR PHYSICIAN.**

Warning Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

- **DO NOT** use this monitor on infants, toddlers, children or persons who cannot express themselves.
- **DO NOT** adjust medication based on readings from this blood pressure monitor. Take medication as prescribed by your physician. **ONLY** a physician is qualified to diagnose and treat high blood pressure and Afib.
- **DO NOT** use this monitor on an injured arm or an arm under medical treatment.
- **DO NOT** apply the arm cuff on your arm while on an intravenous drip or blood transfusion.
- **DO NOT** use this monitor in areas containing high frequency (HF) surgical equipment, magnetic resonance imaging (MRI) equipment, computerized tomography (CT) scanners. This may result in incorrect operation of the monitor and/or cause an inaccurate reading.
- **DO NOT** use this monitor in oxygen rich environments or near flammable gas.
- Consult with your physician before using this monitor if you have common arrhythmias such as atrial or ventricular premature beats or atrial fibrillation; arterial sclerosis; poor perfusion; diabetes; pregnancy; pre-eclampsia or renal disease. **NOTE** that any of these conditions in addition to patient motion, trembling, or shivering may affect the measurement reading.
- **NEVER** diagnose or treat yourself based on your readings. **ALWAYS** consult with your physician.
- To help avoid strangulation, keep the air tube and AC adapter cable away from infants, toddlers and children.
- This product contains small parts that may cause a choking hazard if swallowed by infants, toddlers and children.

Data Transmission

- This product emits radio frequencies (RF) in the 2.4 GHz band. **DO NOT** use this product in locations where RF is restricted, such as on an aircraft or in hospitals. Turn off the Bluetooth® feature in this monitor, remove batteries and/or unplug the AC adapter when in RF restricted areas.

AC Adapter Handling and Usage

- **DO NOT** use the AC adapter if this monitor or the AC adapter cable is damaged. If this monitor or the cable is damaged, turn off the power and unplug the AC adapter immediately.

- Plug the AC adapter into the appropriate voltage outlet. **DO NOT** use in a multi-outlet plug.
- **NEVER** plug in or unplug the AC adapter from the electric outlet with wet hands.
- **DO NOT** disassemble or attempt to repair the AC adapter.

Battery Handling and Usage

- Keep batteries out of the reach of infants, toddlers and children.

Caution

Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient, or cause damage to the equipment or other property.

- Stop using this monitor and consult with your physician if you experience skin irritation or discomfort.
- Consult with your physician before using this monitor on an arm where intravascular access or therapy, or an arterio-venous (A-V) shunt, is present because of temporary interference to blood flow which could result in injury.
- Consult with your physician before using this monitor if you have had a mastectomy.
- Consult with your physician before using this monitor if you have severe blood flow problems or blood disorders as cuff inflation can cause bruising.
- **DO NOT** take measurements more often than necessary because bruising, due to blood flow interference, may occur.
- **ONLY** inflate the arm cuff when it is applied on your upper arm.
- Remove the arm cuff if it does not start deflating during a measurement.
- When the monitor malfunctions, it may get hot. **DO NOT** touch the monitor if it occurs.
- **DO NOT** use this monitor for any purpose other than measuring blood pressure and/or detecting the possibility of Afib.
- During measurement, make sure that no mobile device or any other electrical device that emits electromagnetic fields is within 30 cm of this monitor. This may result in incorrect operation of the monitor and/or cause an inaccurate reading.
- **DO NOT** disassemble or attempt to repair this monitor or other components. This may cause an inaccurate reading.
- **DO NOT** use in a location where there is moisture or a risk of water splashing this monitor. This may damage this monitor.
- **DO NOT** use this monitor in a moving vehicle such as in a car or on an aircraft.
- **DO NOT** drop or subject this monitor to strong shocks or vibrations.
- **DO NOT** use this monitor in places with high or low humidity or high or low temperatures. Refer to section 10.
- During measurement, observe the arm to ensure that the monitor is not causing prolonged impairment to blood circulation.

Important Safety Information

- **DO NOT** use this monitor in high-use environments such as medical clinics or physician offices.
- **DO NOT** use this monitor with other medical electrical (ME) equipment simultaneously. This may result in incorrect operation of the devices and/or cause an inaccurate reading.
- Avoid bathing, drinking alcohol or caffeine, smoking, exercising and eating for at least 30 minutes before taking a measurement.
- Rest for at least 5 minutes before taking a measurement.
- Remove tight-fitting or thick clothing from your arm while taking a measurement.
- Remain still and **DO NOT** talk while taking a measurement.
- **ONLY** use the arm cuff on persons whose arm circumference is within the specified range of the cuff.
- Ensure that this monitor has acclimated to room temperature before taking a measurement. Taking a measurement after an extreme temperature change could lead to an inaccurate reading. OMRON recommends waiting for approximately 2 hours for the monitor to warm up or cool down when the monitor is used in an environment within the temperature specified as operating conditions after it is stored either at the maximum or at the minimum storage temperature. For additional information on operating and storage/transport temperature, refer to section 10.
- **DO NOT** use this monitor after the durable period has ended. Refer to section 10.
- **DO NOT** crease the arm cuff or the air tube excessively.
- **DO NOT** fold or kink the air tube while taking a measurement. This may cause an injury by interrupting blood flow.
- To unplug the air plug, pull on the plastic air plug at the base of the tube, not the tube itself.
- **ONLY** use the AC adapter, arm cuff, batteries and accessories specified for this monitor. Use of unsupported AC adapters, arm cuffs and batteries may damage and/or may be hazardous to this monitor.
- **ONLY** use the approved arm cuff for this monitor. Use of other arm cuffs may result in incorrect readings.
- Inflating to a higher pressure than necessary may result in bruising of the arm where the cuff is applied. **NOTE:** refer to “If your systolic pressure is more than 210 mmHg” in section 3 of this instruction manual for additional information.
- Read and follow the “Correct Disposal of This Product” in section 11 when disposing of the device and any used accessories or optional parts.

Data Transmission

- **DO NOT** replace batteries or unplug the AC adapter while your readings are being transferred to your smart device. This may result in incorrect operation of this monitor and failure to transfer your blood pressure data.

AC Adapter Handling and Usage

- Fully insert the AC adapter into the outlet.
- When unplugging the AC adapter from the outlet, be sure to safely pull from the AC adapter. **DO NOT** pull from the AC adapter cable.
- When handling the AC adapter cable:
 - DO NOT damage it.
 - DO NOT break it.
 - DO NOT tamper with it.
 - DO NOT forcibly bend or pull it.
 - DO NOT twist it.
 - DO NOT use it if it is gathered in a bundle.
 - DO NOT pinch it.
 - DO NOT place it under heavy objects.
- Wipe any dust off of the AC adapter.
- Unplug the AC adapter when not in use.
- Unplug the AC adapter before cleaning this monitor.

Battery Handling and Usage

- **DO NOT** insert batteries with their polarities incorrectly aligned.
 - **ONLY** use 4 “AA” alkaline or manganese batteries with this monitor. **DO NOT** use other types of batteries. **DO NOT** use new and used batteries together. **DO NOT** use different brands of batteries together.
 - Remove batteries if this monitor will not be used for a long period of time.
 - If battery fluid should get in your eyes, immediately rinse with plenty of clean water. Consult with your physician immediately.
 - If battery fluid should get on your skin, wash your skin immediately with plenty of clean, lukewarm water. If irritation, injury or pain persists, consult with your physician.
 - **DO NOT** use batteries after their expiration date.
 - Periodically check batteries to ensure they are in good working condition.
-


1. Know Your Monitor

1.1 Contents/Product Includes

Monitor, arm cuff (HEM-FL31), AC adapter (HHP-CM01), 4 “AA” batteries, storage case, instruction manual, setup instructions, guarantee card, blood pressure pass


1.2 Monitor and Arm Cuff

Front of the monitor

- (a) Display
- (b) [Memory] button
- (c) [Morning/Evening average] button
- (d) [Forward/Backward] buttons
- (e) [User ID selection] switch
- (f) [START/STOP] button


Left and right side of the monitor

g Air jack

h AC adapter jack

Back and bottom of the monitor

i [Connection] button
(Date/Time)

j Battery compartment


1. Know Your Monitor

Arm cuff


(Upper arm circumference: 22 - 42 cm)

Ⓚ Air plug

Ⓛ Air tube

1.3 Display and Symbols

- Ⓐ Evening average symbol
- Ⓑ Morning average symbol
- Ⓒ Memory symbol
- Ⓓ Systolic blood pressure reading
- Ⓔ Comparison display
- Ⓕ Pulse display / Memory number
- Ⓖ Deflation symbol

- Ⓗ Diastolic blood pressure reading
- Ⓘ Average value symbol
- Ⓝ User ID symbol
- Ⓚ Bluetooth ON / OFF symbol
- Ⓛ Sync symbol
- Ⓜ Low battery symbol / Depleted battery symbol
- Ⓝ Date / Time display / Consecutive measurement number for Afib

EN


1. Know Your Monitor

- Ⓒ Afib mode symbol / Afib indicator symbol
- Ⓓ OK symbol
- Ⓔ Hypertension symbol / Morning hypertension symbol
- Ⓕ Movement error symbol
- Ⓖ Irregular heartbeat symbol
- Ⓖ Cuff wrap guide symbol (loose)
- Ⓒ Pulse display / Interval counter for Afib
- Ⓓ Heartbeat symbol
- Ⓖ Cuff wrap guide symbol (OK)

Ⓐ	 AVG	Evening average symbol Appears when viewing the weekly average readings taken in the evening over the past 7 weeks.
Ⓑ	 AVG	Morning average symbol Appears when viewing the weekly average readings taken in the morning over the past 7 weeks.
Ⓒ		Memory symbol Appears when viewing readings stored in the memory.
Ⓓ	Systolic blood pressure reading	
Ⓔ	Comparison display Displays past readings. Refer to section 5 for more information.	
Ⓕ	Pulse display / Memory number Pulse rate appears after the measurement. When pressing the  button, the memory number appears for approximately 1.5 seconds before the pulse rate appears on the display.	
Ⓖ		Deflation symbol Appears during cuff deflation.
Ⓖ	Diastolic blood pressure reading	

①	 AVG	Average value symbol Flashes when viewing the average of the latest 2 or 3 readings taken within the most recent 10 minute timeframe.
②		User ID symbol The user ID number appears when selecting the user with the [user ID selection] switch.
③		Bluetooth ON symbol Appears when your readings are being transferred.
		Bluetooth OFF symbol Appears when Bluetooth is disabled.
④		Sync symbol Flashes/appears when your data needs to be transferred because the stored memory is either almost, or completely full. Once you pair your monitor with your smart device, transfer your readings immediately before the monitor deletes the oldest reading. Up to 100 readings can be stored in the internal memory of your monitor per user.
⑤		Low battery symbol Flashes when batteries are low.
		Depleted battery symbol Appears when batteries are depleted.
⑥	Date / Time display / Consecutive measurement number for Afib	

1. Know Your Monitor

		<p>Afib mode symbol Appears when a measurement is taken if the Afib mode is enabled. The Afib mode automatically takes 3 consecutive readings at 30-second intervals and displays the average. To use the Afib mode, refer to sub-section 3.1. To view your readings in Afib measurement, refer to sub-section 5.1.</p>
Ⓞ		<p>Afib indicator symbol Appears if a possibility of atrial fibrillation (Afib) was detected during the Afib mode measurement. In the Afib mode, this particular irregularity of the heart rhythm is indicated if it deviates +/- 25 % than the average rhythm detected while your monitor is measuring blood pressure. If it continues to appear, we recommend you to consult with and follow the directions of your physician.</p>
Ⓟ		<p>OK symbol Flashes when your monitor is connected to your smart device or readings are transferred successfully.</p>
Ⓞ	 	<p>Hypertension symbol Appears if your systolic blood pressure is 135 mmHg or above and/or the diastolic blood pressure is 85 mmHg or above.</p> <p>Morning hypertension symbol Appears when the morning weekly average is 135/85 mmHg or above.</p>

<p>Ⓡ</p>		<p>Movement error symbol Appears with your reading when your body is moving during a measurement. The body movement function is disabled when the possibility of Afib or irregular heartbeat is detected during a Afib mode measurement. If it appears, remove the arm cuff and wait for 2 to 3 minutes. Then, reapply the cuff, remain still and take another measurement.</p>
<p>Ⓢ</p>		<p>Irregular heartbeat symbol Appears when irregular rhythm is detected 2 or more times during a measurement. It does not appear during a Afib mode measurement.</p> <p>An irregular heartbeat rhythm is defined as a rhythm that is 25 % less or 25 % more than the average rhythm detected while your monitor is measuring blood pressure.</p> <div style="text-align: center;">  <p>Irregular Heartbeat Short Long</p> <p>Normal Heartbeat</p> </div> <p>Irregular heartbeat symbol may be displayed when the monitor detects arm movements. Refer to section 7 for more information.</p> <p>If it continues to appear, we recommend you to consult with and follow the directions of your physician.</p>
<p>Ⓣ</p>		<p>Cuff wrap guide symbol (loose) Appears if the arm cuff is not wrapped around the arm correctly while taking a measurement. It also appears when viewing past readings.</p>
<p>Ⓤ</p>		<p>Pulse display / Interval counter for Afib Pulse rate appears after the measurement. Counts the seconds for next Afib measurement during Afib mode.</p>

1. Know Your Monitor

		Heartbeat symbol Flashes while taking a measurement.
		Cuff wrap guide symbol (OK) Appears if the arm cuff is wrapped around the arm correctly while taking a measurement. It also appears when viewing past readings.

The difference between Afib indicator function and irregular heartbeat function

Afib indicator function: Detects Afib possibility in 3-times measurement. (based on the clinical paper*)

Irregular heartbeat function: Detects irregular heartbeat including Afib in 1 measurement.

* M. Ishizawa, T. Noma, T. Minamino et al. Multiple measurements with automated blood pressure monitor can detect atrial fibrillation with high sensitivity and specificity in general cardiac patients, ESC Congress 2018

2018 ESH/ESC** Guidelines for the management of arterial hypertension

Definitions of hypertension by office and home blood pressure levels

	Office	Home
Systolic Blood Pressure	≥ 140 mmHg	≥ 135 mmHg
Diastolic Blood Pressure	≥ 90 mmHg	≥ 85 mmHg

These ranges are from statistical values for blood pressure.

** European Society of Hypertension (ESH) and European Society of Cardiology (ESC).


⚠ Warning

• **NEVER** diagnose or treat yourself based on your readings. **ALWAYS** consult with your physician.

2. Preparation

2.1 Installing Batteries

1. Push down on the hook of the battery cover and pull upward.

EN

2. Insert 4 “AA” batteries as indicated in the battery compartment.
3. Close the battery cover.

As soon as batteries are correctly installed, the year flashes on the display.

Refer to sub-section 2.3 to set the date and time manually, or refer to sub-section 4.1 to pair your monitor to set the date and time automatically.

Note

- When the “” symbol flashes on the display, it is recommended you replace batteries.
- To replace batteries, turn your monitor off and remove all batteries. Then replace with 4 new batteries at the same time.
- Replacing batteries will not delete previous readings.
- The supplied batteries may have a shorter life span than new batteries.
- Disposal of used batteries should be carried out in accordance with local regulations.


EN18

2.2 Using the AC Adapter

Note

- Make sure not to place your monitor in a location where it is difficult to plug and unplug the AC adapter.
- We recommend keeping batteries in your monitor at all times, even if you choose to use the AC adapter. If only the AC adapter is used without keeping batteries in your monitor, you may need to reset the date and time each time you unplug and plug back the AC adapter. The readings will not be deleted.

- 1. Insert the AC adapter plug into the AC adapter jack on the right side of your monitor.**

- 2. Plug the AC adapter into an electrical outlet.**

To unplug the AC adapter, unplug the AC adapter from the electrical outlet and then remove the AC adapter plug from the monitor.

⚠ Caution

- **ONLY use the AC adapter, arm cuff, batteries and accessories specified for this monitor. Use of unsupported AC adapters, arm cuffs and batteries may damage and/or may be hazardous to this monitor.**

2.3 Setting the Date and Time


Note

- Set your monitor to the correct date and time before taking a measurement for the first time.
- The date and time will automatically be set when the monitor is paired with “OMRON connect” app. Refer to sub-section 4.1.
- As soon as batteries are correctly installed, the year flashes on the display of your monitor.

EN

1. If nothing is displayed, press the button (rear surface of the monitor) to display the year.

If the year is already displayed, skip this step.


2. Press the or button to adjust the year.

Press and hold the  button to quickly advance the years.

Press and hold the  button to quickly go backwards through the years.

3. Press the  button to confirm the year, then the month will flash. Repeat the same steps to adjust the month, day, hour, and minute.

4. Press the  button to store the setting.

Your monitor automatically turns off 3 seconds after the minutes have been set in this sequence.

Note

- If you need to reset the date and time, replace batteries or press the  button, then set the date and time as shown in above steps.
- If the date and time are not set, “-:--” flashes during or after measurement.

2.4 Blood Pressure Measurement Tips

To help ensure an accurate reading, follow these directions:

- Stress raises blood pressure. Avoid taking measurements during stressful times.
- Measurements should be taken in a quiet place.
- It is important to take measurements at the same times each day. Taking measurements in the morning and in the evening is recommended.
- Remember to have a record of your blood pressure and pulse readings for your physician. A single measurement does not provide an accurate indication of your true blood pressure. Please use the included Blood Pressure Pass to keep records of several readings over a certain period of time. You may also download the PDF files of the pass by visiting www.omron-healthcare.com.

EN

Caution

- **Avoid bathing, drinking alcohol or caffeine, smoking, exercising and eating for at least 30 minutes before taking a measurement.**
 - **Rest for at least 5 minutes before taking a measurement.**
-

2.5 Applying the Arm Cuff


Note


- The following steps are for applying the arm cuff to your left arm. When you take a measurement on your right arm, follow the cuff wrap instructions for a right arm use mentioned under “To take a measurement on your right arm” located at the end of this sub-section.
- The blood pressure can differ between the right arm and the left arm, and the measured blood pressure readings can be different. OMRON recommends to always use the same arm for measurement. If the readings between both arms differ substantially, check with your physician to determine as to which arm to use for your measurements.

Caution

- Remove tight-fitting or thick clothing from your arm while taking a measurement.

1. **Plug the arm cuff into your monitor by inserting the air plug into the air jack securely until it clicks.**


2. Apply the arm cuff to your left upper arm.

EN

Note

- The bottom edge of the arm cuff should be 1 to 2 cm above the inside elbow. The air tube is on the inside of your arm and aligned with your middle finger.

3. Make sure that the air tube is positioned on the inside of your arm and wrap the cuff securely, so it can not move around your arm.**To take a measurement on your right arm**

When you take a measurement on your right arm, the air tube should run along the side of your elbow, along the bottom of your arm. Be careful not to rest your arm on the air tube.

2.6 Sitting Correctly

To take a measurement, you need to be relaxed and comfortably seated in a room with a comfortable temperature. Place your arm on the table.


- Sit in a comfortable chair with your back and arm supported.
- Keep your feet flat and your legs uncrossed.
- The arm cuff should be placed on your arm at the same level as your heart, with the arm resting comfortably on a table.

2.7 Selecting a User ID

Your monitor is designed to take measurements and store the readings in the memory for 2 people using the [user ID selection] switch.

1. Select your user ID (1 or 2).

3. Using Your Monitor without a Smart Device

3.1 Taking a Measurement

You can use your monitor without pairing to a smart device. To pair your monitor with a smart device, refer to sub-section 4.1.

Note

- To stop the measurement, press the [START/STOP] button once to deflate the arm cuff.

⚠ Caution

- **DO NOT use this monitor with other medical electrical (ME) equipment simultaneously. This may result in incorrect operation of the devices and/or cause an inaccurate reading.**
- **Remain still and DO NOT talk while taking a measurement.**

EN

1. Press the [START/STOP] button.

All symbols appear on the display before starting the measurement.

2. Remain still and do not talk until the entire measurement process is complete.

The “” symbol flashes at every heartbeat.


Note

- The “” symbol appears if the arm cuff is wrapped around the arm correctly.
- When the “” symbol is displayed, the arm cuff is not applied correctly. Press the [START/STOP] button to turn your monitor off, then apply the cuff correctly.

3. Using Your Monitor without a Smart Device

After your monitor has detected your blood pressure and pulse rate, the cuff automatically deflates. Your blood pressure and pulse rate are displayed with the previous reading.

If either the systolic or the diastolic reading is high (refer to sub-section 1.3), the “” symbol appears.


3. Press the [START/STOP] button to turn your monitor off.

Note

- Your monitor automatically turns off after 2 minutes.
- Wait 2 - 3 minutes between measurements. The wait time allows the arteries to decompress and return to their pre-measurement state. You may need to increase the wait time depending on your individual physiological characteristics.

■ Using the Afib Mode

Afib indicator function

In the Afib mode, the monitor automatically takes 3 consecutive readings at 30-second intervals and displays the average. If there is a possibility of Afib, the “” symbol will appear.

This is not the diagnosis, it is only a potential finding for Afib.

You should contact your physician to discuss the findings.

If you are experiencing any symptoms, contact a medical professional.

EN


1. **Select your user ID.**
2. **Press and hold the [START/STOP] button for more than 2 seconds.**
3. **Release the [START/STOP] button when “AF” appears on the display.**

The arm cuff starts to inflate automatically.

After the first measurement is complete, wait for the next measurement to start.


The 30-second interval counter will appear on the display.


3. Using Your Monitor without a Smart Device

After the second measurement is complete, wait for the next measurement to start.


The 30-second interval counter will appear on the display.

After the third measurement is complete, the average for the 3 measurements appears on the display. If Afib is detected, the “” symbol will appear.

Note

- If error symbol “” appears, follow the instructions in section 7.

To view readings in Afib measurement, refer to sub-section 5.1.

■ Using the Guest Mode

Your monitor stores readings for 2 users in the memory. The guest mode can be used to take a single measurement for another user.

Note

- No readings are stored in the memory, and the Afib mode is not available when the guest mode is selected.

1. While holding the  button down, press the [START/STOP] button.

The user ID symbol and the date/time display will not be shown.

EN

2. Release the  and [START/STOP] buttons.

The arm cuff automatically starts to inflate.

■ If your systolic pressure is more than 210 mmHg

If your systolic pressure is more than 210 mmHg, take a measurement manually.

After the arm cuff starts to inflate, press and hold the [START/STOP] button until the monitor inflates 30 to 40 mmHg higher than your expected systolic pressure.

Note

- Do not inflate above 299 mmHg.

Caution

- Inflating to a higher pressure than necessary may result in bruising of the arm where the cuff is applied.

4. Using Your Monitor with a Smart Device

4.1 Pairing Your Monitor with a Smart Device

The date and time on your monitor will automatically be set when you pair your monitor with your smart device.

Review the list of compatible smart devices at www.omronconnect.com/devices/


1. **Enable Bluetooth on your smart device.**
2. **Download and install the free “OMRON connect” app onto your smart device.**

If you already have the “OMRON connect” app and have created your account, open the app and add your new monitor.

3. **Open the app and follow the pairing instructions shown on your smart device.**
4. **Confirm that your monitor is connected successfully.**


When your monitor is connected successfully to your smart device, the OK symbol appears.

5. Press the [START/STOP] button to turn your monitor off.

Note

- If the “Err” symbol appears, follow the instructions in the “OMRON connect” app.

EN

- Your monitor automatically turns off after 2 minutes.
- Be aware that OMRON will not be responsible for the loss of data and/or information in the app.
- “OMRON connect” is the only app that we recommend to use with your monitor to transfer data correctly.

4.2 Taking a Measurement

Refer to sub-section 3.1 to start a measurement. Open the app on your smart device and follow the instructions.

Note

- If your profile has been registered, once the [START/STOP] button is pressed, your birthday is displayed and the user ID symbol flashes twice.
- Profile registration is available on the “OMRON connect” app. Refer to sub-section 4.1.
- Make sure that your monitor is within 5 m of your smart device.

■ Transferring Data Automatically/Manually

- When Bluetooth is enabled, the readings will be transferred to your smart device within 1 hour after taking a measurement.
- To transfer the readings manually, press the  button.

4. Using Your Monitor with a Smart Device

4.3 Disabling/Enabling Bluetooth


Disable Bluetooth in your monitor in the following areas where use of wireless equipment is prohibited.

- On aircraft
- In hospitals
- While abroad

Bluetooth is enabled as default.

1. **When your monitor is off, press and hold the  button (rear surface of the monitor) for more than 10 seconds.**


“oFF” flashes 3 times.

2. **Press the [START/STOP] button to turn your monitor off.**

Note

- To enable Bluetooth, press and hold the  button for more than 2 seconds. “on” appears on the display.
- Your monitor automatically turns off 2 minutes after the completion.
- When Bluetooth is disabled, the “ off” symbol appears.

5. Using Memory Function

Your monitor automatically stores up to 100 readings per user.

Note

- This monitor can be used by multiple users. Refer to sub-section 2.7 to switch user. Make sure which user ID stores your own readings.


5.1 Viewing the Readings Stored in Memory


Your monitor has “comparison mode” view function, that shows you the latest reading and the past reading at the same time.

1. Select your user ID.

2. Press the button.

The readings appear in comparison mode as below.

If the past reading was taken in Afib mode, the “3” and “” symbols or the only “3” symbol appear. The reading shown on the left is the average of 3 measurements.


If the latest reading was taken in Afib mode, the “3” and “” symbols or the only “3” symbol appear. The reading shown on the right is the average of 3 measurements.

After a memory number appears, the latest reading will appear with a past reading. A past reading will appear in the comparison display on the left.

5. Using Memory Function

Note

- If the reading is high (refer to sub-section 1.3), the “” symbol appears.
- If there are no readings stored in the memory, the screen to the right is displayed.
- If the memory is full, new readings will replace the old readings, starting with the oldest reading.

3. Press the or button repeatedly to scroll through the previous readings stored in the memory.

 : To view the older readings

 : To view the more recent readings

5.2 Viewing the Average Reading

Your monitor can calculate an average reading based on the latest 2 or 3 readings taken within the most recent 10 minute timeframe.

Note

- Set your monitor to the correct date and time before taking a measurement to view the average readings.


1. Select your user ID.

2. When your monitor is off, press and hold the button for more than 3 seconds.

The “ AVG ” symbol flashes.

Note

- If there are only 2 readings in the memory for that period, the average will be based on these 2 readings.
- If the average reading is high (refer to sub-section 1.3), the “” symbol appears.

5.3 Viewing Weekly Averages

Your monitor calculates and displays weekly averages for the measurements taken in the morning and evening for 7 weeks per user.

Note


- Set your monitor to the correct date and time before taking a measurement to view the weekly average readings.

Morning Weekly Average Calculation

This is the average for the measurements taken during the morning (4:00 - 9:59) between Sunday and Saturday. The 2 or 3 readings taken within the first 10 minute timeframe in the morning between 4:00 - 9:59 will be used to calculate the morning average for each day.

Evening Weekly Average Calculation

This is the average for the measurements taken during the evening (19:00 - 1:59) between Sunday and Saturday. The 2 or 3 readings taken within the last 10 minute timeframe in the evening between 19:00 - 1:59 will be used to calculate the evening average for each day.



Note


- When the measurements taken from Thursday, the average is calculated between Thursday and Saturday.

1. Select your user ID.

2. Press the button once.

For the morning weekly average for the current week, the “ AVG” symbol and “THIS WEEK” appear on the display. A past reading will appear in the comparison display on the left. If the morning weekly average is high (refer to sub-section 1.3), the “” or “” symbol appears.


Press the  button again for the evening weekly average.

The “ AVG” symbol and “THIS WEEK” appear on the display. A past reading will appear in the comparison display on the left.

If the evening weekly average is high, the “” or “” symbol appears.

3. Press the or button repeatedly to scroll through the previous readings stored in the memory.

Your monitor displays “1 wk[s]” for the previous week's average, to “7 wk[s]” for the oldest set of averages.

4. Press the [START/STOP] button to turn off your monitor.

6. Other Settings

6.1 Deleting All Readings Stored in Memory


1. Select your user ID.

2. Press the  button.

After a memory number appears, the latest reading will appear with a past reading.


3. While holding the  button down, press and hold the [START/STOP] button for more than 3 seconds.

The display as below flashes.

4. While holding the  button down again, press and hold the [START/STOP] button for more than 3 seconds.

The display stops flashing, then all readings are deleted.

EN

5. Press the [START/STOP] button to turn your monitor off.

Note


- All readings will be deleted. You cannot partially delete the readings stored in the memory.
- Your monitor automatically turns off after 2 minutes.

6.2 Restoring Your Monitor to the Default Settings

To delete all the information stored in your monitor, follow the instructions below. Make sure that your monitor is turned off.


1. While holding the button (rear surface of the monitor) down, press and hold the [START/STOP] button for more than 3 seconds.

The display as below flashes.

2. While holding the button down again, press and hold the [START/STOP] button for more than 3 seconds.

The display stops flashing, then your monitor is restored to the default settings.

3. Press the [START/STOP] button to turn your monitor off.

Note

- After the completion, the year flashes on the display. Refer to sub-section 2.3 to set the date and time manually, or refer to sub-section 4.1 to pair your monitor to set the date and time automatically.
- Reverting to the default setting of your monitor does not delete the information in the app.

EN


7. Error Messages and Troubleshooting

If any of the below problems occur during measurement, check to make sure that no other electrical device is within 30 cm. If the problem persists, refer to the table below.

Display/ Problem	Possible Cause	Solution
E1 appears or the arm cuff does not inflate.	The [START/STOP] button was pressed while the arm cuff is not applied.	Press the [START/STOP] button again to turn the monitor off. After inserting the air plug securely and applying the arm cuff correctly, press the [START/STOP] button.
	Air plug is not completely plugged into the monitor.	Insert the air plug securely. Refer to sub-section 2.5.
	The arm cuff is not applied correctly.	Apply the arm cuff correctly, then take another measurement. Refer to sub-section 2.5.
	Air is leaking from the arm cuff.	Replace the arm cuff with a new one. Contact your OMRON retail outlet or distributor.

Display/ Problem	Possible Cause	Solution
<p>E2 appears or a measurement cannot be complete after the arm cuff inflates.</p>	<p>You move or talk during a measurement and the arm cuff does not inflate sufficiently.</p>	<p>Remain still and do not talk during a measurement. If “E2” appears repeatedly, inflate the arm cuff manually until the systolic pressure is 30 to 40 mmHg above your previous readings. Refer to sub-section 3.1.</p>
	<p>The systolic pressure is above 210 mmHg and a measurement cannot be taken.</p>	<p>Inflate the arm cuff manually until the systolic pressure is 30 to 40 mmHg above your previous readings. Refer to sub-section 3.1.</p>
<p>E3 appears</p>	<p>The arm cuff is inflated exceeding the maximum allowable pressure.</p>	<p>Do not touch the arm cuff and/or bend the air tube while taking a measurement. If inflating the arm cuff manually, refer to the end of sub-section 3.1.</p>
<p>E4 appears</p>	<p>You move or talk during a measurement. Vibrations disrupt a measurement.</p>	<p>Remain still and do not talk during a measurement.</p>

7. Error Messages and Troubleshooting

Display/ Problem	Possible Cause	Solution
 appears	The pulse rate is not detected correctly.	Apply the arm cuff correctly, then take another measurement. Refer to sub-section 2.5. Remain still and sit correctly during a measurement. If the “  ” or “  ” symbol continues to appear, we recommend you to consult with your physician.
 appears		
 appears		
 appears		
 does not flash during a measurement		
 appears	Blood pressure measurements were not taken correctly in an Afib mode measurement.	
 appears	The monitor has malfunctioned.	Press the [START/STOP] button again. If “Er” still appears, contact your OMRON retail outlet or distributor.

Display/ Problem	Possible Cause	Solution
 appears	The monitor cannot connect to a smart device or transmit data correctly.	Follow the instructions shown in the “OMRON connect” app. If the “Err” symbol still appears after checking the app, contact your OMRON retail outlet or distributor.
 flashes	The monitor is waiting for pairing with the smart device.	Refer to sub-section 4.1 for pairing your monitor with your smart device, or press the [START/STOP] button to cancel pairing and turn your monitor off.
 flashes	The monitor is ready to transfer your readings to the smart device.	Open the “OMRON connect” app to transfer your readings.
 flashes	There are more than 80 readings in memory to be transferred. The date and time is not set.	Pair or transfer your readings to the “OMRON connect” app so you can keep them in memory in the app, and this error symbol disappears.
 appears	There are 100 readings in memory to be transferred.	
 flashes	Batteries are low.	Replacing all 4 batteries with new ones is recommended. Refer to sub-section 2.1.

7. Error Messages and Troubleshooting

Display/ Problem	Possible Cause	Solution
 appears or the monitor is turned off unexpectedly during a measurement.	Batteries are depleted.	Immediately replace all 4 batteries with new ones. Refer to sub-section 2.1.
Nothing appears on the display of the monitor.	Batteries are depleted.	Check the battery installation for proper placement. Refer to sub-section 2.1.
	Battery polarities are not properly aligned.	
Readings appear too high or too low.	Blood pressure varies constantly. Many factors including stress, time of day, and/or how you apply the arm cuff, may affect your blood pressure. Review sub-sections 2.4 - 2.6 and section 3.	
Any other communication issue occurs.	Follow the instructions shown in the smart device, or visit the “Help” section in the “OMRON connect” app for further help. If the problem still persists, contact your OMRON retail outlet or distributor.	
Any other problems occur.	Press the [START/STOP] button to turn the monitor off, then press it again to take a measurement. If the problem continues, remove all batteries and wait for 30 seconds. Then re-install batteries. If the problem still persists, contact your OMRON retail outlet or distributor.	

Frequently Asked Questions regarding Afib

<p>What is different between the Afib indicator function and ECG?</p>	<p>The Afib indicator function and ECG use different technologies. An ECG measures the electrical activity of the heart and can be used to diagnose Afib. The Afib indicator function detects irregular heartbeat and can suggest the possibility of Afib with a sensitivity of 95.5 % and specificity of 93.8 %. Refer to section 12 for details.</p>
<p>If the “” symbol does not appear, it means there is no possibility of Afib?</p>	<p>Even if the “” symbol does not appear, there is still a possibility of Afib.</p>
<p>Should I consult with my physician if I the “” symbol appears?</p>	<p>We recommend you to consult with your physician because there is a possibility of Afib. However, the “” symbol may be displayed for other reasons, such as other heart arrhythmias.</p>
<p>What is different between Afib indicator function and irregular heartbeat function?</p>	<p>The irregular heartbeat function detects irregularities in the pulse waves in one measurement. The Afib indicator function suggests the possibility of Afib when blood pressure is measured 3 consecutive times.</p>
<p>What should I do if the “” symbol sometimes appears?</p>	<p>Afib does not always have symptoms. We recommend you to consult with and follow the directions of your physician.</p>

7. Error Messages and Troubleshooting

<p>I have been diagnosed with Afib by the physician, but the “” symbol does not appear.</p>	<p>Afib may not occur at the time of specific blood pressure measurements. We recommend you to consult with your physician regularly.</p>
<p>Is the blood pressure reading reliable when the “” symbol appears?</p>	<p>Afib or an irregular heartbeat can influence your blood pressure measurements and make it difficult to get an accurate reading. Repeated measurements may be required to overcome variabilities.* In Afib mode, the blood pressure measurement is taken 3 times, and the average is displayed. The monitor will indicate an error message (E5/E6) if the influence of the irregular heartbeat is too severe to give a measurement result. If this occurs repeatedly, we recommend that you consult with your physician.</p> <hr/> <p>* Prof. Roland Asmar et al. European Society of Hypertension Recommendations for Conventional, Ambulatory and Home Blood Pressure Measurement</p>

8. Maintenance

8.1 Maintenance

To protect your monitor from damage, follow the directions below:

- Changes or modifications not approved by the manufacturer will void the user warranty.
- The monitor does not require special maintenance.

Caution

- **DO NOT** disassemble or attempt to repair this monitor or other components. This may cause an inaccurate reading.

EN

8.2 Storage

Keep your monitor and other components in the storage case when not in use.


Store your monitor and other components in a clean, safe location.

1. Remove the arm cuff from the monitor.

Caution

- To unplug the air plug, pull on the plastic air plug at the base of the tube, not the tube itself.

2. Gently fold the air tube into the arm cuff.

Note

- Do not bend or crease the air tube excessively.

3. Place your monitor and other components in the storage case.

Do not store your monitor and other components:

- If your monitor and other components are wet.
- In locations exposed to extreme temperatures, humidity, direct sunlight, dust or corrosive vapors such as bleach.
- In locations exposed to vibrations or shocks.

To protect your monitor during storage, an optional LCD cover is available as accessory. Refer to section 9.

8.3 Cleaning

- Do not use any abrasive or volatile cleaners.
- Use a soft dry cloth or a soft cloth moistened with mild (neutral) detergent to clean your monitor and arm cuff and then wipe them with a dry cloth.
- Do not wash or immerse your monitor and arm cuff or other components in water.
- Do not use gasoline, thinners or similar solvents to clean your monitor and arm cuff or other components.


8.4 Calibration and Service

- The accuracy of this blood pressure monitor has been carefully tested and is designed for a long service life.
- It is generally recommended to have the unit inspected every two years to ensure correct functioning and accuracy. Please consult your authorised OMRON dealer or the OMRON Customer Service at the address given on the packaging or attached literature.

9. Description of accessories, medical devices, or non-medical devices intended for use in combination with the automatic blood pressure monitor OMRON M7 Intelli IT.


Arm cuff

Arm circumference
22 - 42 cm

HEM-FL31

AC Adapter


HHP-CM01

EN

Note

- Do not throw the air plug away. The air plug can be applicable to the optional cuff.

Protective LCD Cover

Model: HEM-CACO-734

For further information, visit our website at www.csmedica.ru

10. Specifications

Product description	OMRON Automatic Upper Arm Blood Pressure Monitor
Model (code)	M7 Intelli IT (ALRU)
Display	LCD Digital Display
Measurement method	Oscillometric method
Cuff pressure range	From 0 to 299 mmHg
Blood pressure measurement range	From 40 to 260 mmHg***
Pulse measurement range	From 40 to 180 beats / min.
Accuracy	Pressure: ± 3 mmHg Pulse: ± 5 % of display reading
Inflation	Automatic by electric pump
Deflation	Automatic pressure release valve
Transmission method	Bluetooth® Low Energy
Wireless communication	Frequency range: 2.4 GHz (2400 - 2483.5 MHz) Modulation: GFSK Effective radiated power: < 20 dBm
Operation mode	Continuous operation
Memory	100 measurements per user
Rating	DC6 V 4.0 W
Power source	4 "AA" batteries 1.5 V or AC adapter (INPUT AC100 - 240 V 50 - 60 Hz 0.12 - 0.065 A)
Durable period	Monitor: 5 years Arm cuff: 5 years AC adapter: 5 years
Battery life	Approx. 1000 measurements (using new alkaline batteries) The number of times may decrease when using Afib mode because one Afib indication consists of 3 regular measurements.

Applied part	Type BF (cuff)
Protection against electric shock	Internally powered ME equipment (When using only batteries) Class II ME equipment (AC adapter)
IP classification	
Monitor - IP classification is degrees of protection provided by enclosures in accordance with IEC 60529. This unit is protected against solid foreign objects of diameter 12.5 mm such as a finger and greater.	IP20
AC adapter - IP classification is degrees of protection provided by enclosures in accordance with IEC 60529. This unit is protected against solid foreign objects of diameter 12.5 mm such as a finger and greater. This unit is protected against the ingress of vertically falling water drips.	IP21
Operating conditions	
Temperature/	From +10 to +40 °C
Humidity/	From 15 to 90 % RH (not condensed)
Air pressure	From 800 to 1060 hPa

10. Specifications

Storage and Transport conditions	
Temperature/ Humidity	From -20 to +60 °C From 10 to 90 % RH (not condensed)
Weight	Monitor: 460 ±46 g without batteries Arm cuff: 161 ±17 g AC adapter: 47.5 ±5 g Storage case: 120 ±12 g
Outer dimensions	Monitor: 191 ±10 x 85 ±5 x 120 ±6 mm (w x h x l) Arm cuff: 532 ±5 x 15 ±3 x 171 ±1 mm (w x h x l) Arm cuff air tube length: 750 ±38 mm AC adapter: 21 ±2 x 54.2 ±3 x 64.5 ±4 mm (w x h x l) AC adapter cord length: 1500 ±75 mm Storage case: 210 ±11 x 108 ±6 x 150 ±8 mm (w x h x l)
Cuff circumference applicable to the monitor	From 22 to 42 cm
Cuff / Tube material	Nylon, polyester, polyvinyl chloride
Contents	Monitor, arm cuff (HEM-FL31), AC adapter (HHP-CM01), 4 "AA" batteries, storage case, instruction manual, setup instructions, guarantee card, blood pressure pass
Maximum temperature of the applied part	Lower than +43 °C

Note

- This device is clinically investigated according to the requirements of ISO 81060-2:2013 (excluding pregnant and pre-eclampsia patients). In the clinical validation study, K5 was used on 85 subjects for determination of diastolic blood pressure.
- This device has been validated for use on pregnant and pre-eclampsia patients according to the Modified European Society of Hypertension Protocol*
- This device has been validated for use on diabetic (Type II) population**

* Topouchian J et al. *Vascular Health and Risk Management* 2018:14 189–197

** Chahine M.N. et al. *Medical Devices: Evidence and Research* 2018:11 11–20

*** SYS: 60 to 260 mmHg
DIA: 40 to 215 mmHg

About a wireless communication interference

This product operates in an unlicensed ISM band at 2.4 GHz. In the event this product is used near other wireless devices such as microwave and wireless LAN, which operate on the same frequency band as this product, there is a possibility that interference may occur.

If interference occurs, stop the operation of the other devices or relocate this product away from other wireless devices before attempting to use it.

11. Guidance and Manufacturer's Declaration

CE 0197

- This blood pressure monitor is designed according to the European Standard EN1060, Non-invasive sphygmomanometers Part 1: General Requirements and Part 3: Supplementary requirements for electromechanical blood pressure measuring systems.
- Hereby, OMRON HEALTHCARE Co., Ltd., declares that the radio equipment type M7 Intelli IT (ALRU) is in compliance with Directive 2014/53/EU.
- The full text of the EU declaration of conformity is available at the following internet address: www.omron-healthcare.com
- This OMRON product is produced under the strict quality system of OMRON HEALTHCARE Co., Ltd., Japan. The Core component for OMRON blood pressure monitors, which is the Pressure Sensor, is produced in Japan.
- Notify the authorized representative of the manufacturer of any serious incidents involving this product.


Applied Standards

EN ISO 15223-1:2016	EN ISO 10993-5:2009
EN 1041:2008+A1:2013	EN ISO 10993-10:2013
EN 1060-1:1995+A2:2009	EN ISO 13485:2016
EN 1060-3:1997+A2:2009	EN ISO 14971:2012
EN 60601-1:2006+A1:2013	EN ISO 81060-2:2014
EN 60601-1-2:2015	EN 300 328 V2.2.2
EN 60601-1-6:2010+A1:2015	EN 301 489-1 V2.2.3
EN 60601-1-11:2015	EN 301 489-17 V3.2.4
EN 80601-2-30:2010+A1:2015	EN 62479:2010
EN 62304:2006+A1:2015	EN IEC 62368-1:2020+A11:2020
EN 62366-1:2015	EN IEC 63000:2018
EN ISO 10993-1:2009/AC:2010	

Symbols description	
	Applied part - Type BF Degree of protection against electric shock (leakage current)
	Class II equipment. Protection against electric shock
IP XX	Ingress protection degree provided by IEC 60529
	CE Marking
	Metrology symbol
	Symbol of Eurasian Conformity
	Serial number
	LOT number
	Reference catalogue number
	Temperature limitation
	Humidity limitation
	Atmospheric pressure limitation
	Indication of connector polarity

EN

11. Guidance and Manufacturer's Declaration

Symbols description	
	For indoor use only
	OMRON's trademarked technology for blood pressure measurement
	Identifier of cuffs compatible for the device
	Artery Mark
	Manufacturer's quality control mark
LATEX FREE	Not made with natural rubber latex
	Arm circumference
	Need for the user to consult this instruction manual.
	Need for the user to follow this instruction manual thoroughly for your safety.
	Direct current
	Alternating current
	To indicate generally elevated, potentially hazardous, levels of non-ionizing radiation, or to indicate equipment or systems e.g. in the medical electrical area that include RF transmitters or that intentionally apply RF electromagnetic energy for diagnosis or treatment.

Symbols description	
	OMRON connect - App for uploading measurement data to smart device.
 ТЕХНОЛОГИЯ OMRON HEALTHCARE ЯПОНИЯ	Technology by OMRON Healthcare in Japan
Product production date is integrated in the Serial number, which is placed on the product and/or sales package: the first 4 digits mean year of production, the next 2 digits mean month of production.	

EN

Correct Disposal of This Product (Waste Electrical & Electronic Equipment)

This marking shown on the product or its literature, indicates that it should not be disposed of, with other household wastes at the end of its working life. To prevent possible harm to the environment or human health from uncontrolled waste disposal, please separate this product from other types of wastes and recycle it responsibly to promote the sustainable reuse of material resources.

Household users should contact either the retailer where they purchased this product, or their local government office, for details of where and how they can return this item for environmentally safe recycling.

Business users should contact their supplier and check the terms and conditions of the purchase contract. This product should not be mixed with other commercial waste for disposal.

DISPOSAL

At the end of the life of the product, it must be disposed in accordance with national and regional regulations. In medical establishments, in accordance with the requirements of SanPin 2.1.3684-21 "Sanitary and epidemiological requirements for the maintenance of the territories of urban and rural settlements, water bodies, potable water and potable water supply, air, soils, residential quarters, operation of industrial, public premises, organization and implementation of sanitary and anti-epidemic (preventive) measures", the product should be recycled as a Class A product in accordance with the requirements of the Sanitary Rules applicable to municipal solid waste management.

Battery Disposal

For battery disposal, please contact special localized collection points located in your city, or local authorities for detailed information on where and how to return the batteries for environmentally safe recycling.

VERIFICATION

The document for verification is P 1323565.2.001-2018.

Reverification interval 2 years.

The type approval mark is labelled on the accompanying documents. In accordance with Article 13(4) of Federal Law No. 102-ФЗ [102-FZ] "On ensuring the uniformity of measurements", the results of the calibration are confirmed through the inclusion of information on the calibration results in the Federal Information Fund for Ensuring the Uniformity of Measurements. Confirmation of the calibration referring to the name, model and serial number of the device can also be found on this website: www.csmedica.ru.

Electromagnetic Compatibility

Important information regarding Electromagnetic Compatibility (EMC)

M7 Intelli IT (ALRU) conforms to EN60601-1-2:2015

Electromagnetic Compatibility (EMC) standard.

Nevertheless, special precautions need to be observed:

- The use of accessories and cables other than those specified or provided by OMRON could result in increased electromagnetic emission or decreased electromagnetic immunity of the device and result in improper operation.
- During measurement, the use of the device adjacent to or stacked with other device should be avoided because it could result in improper operation. In case such use is necessary, the device and other device should be observed to verify that they are operating normally.
- During measurement, Portable RF communications device (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm to any part of the device, including cables specified by OMRON. Otherwise, degradation of the performance of the device could result.

Table 1 - EMISSION Limits and Compliance

Phenomenon	EMISSION Limits	Conformity
Conducted and radiated RF EMISSIONS	CISPR 11	Group1, Class B
Voltage fluctuations and flicker	See IEC 61000-3-3	Complies

Table 2 - IMMUNITY TEST LEVELS

Phenomenon	Basic EMC standard	IMMUNITY TEST LEVELS
Electrostatic discharge	IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air for enclosure port
Radiated RF electromagnetic fields	IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz 80 % AM at 1 kHz for enclosure port
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	See table 3
Electrical fast transients / bursts	IEC 61000-4-4	±2 kV for Input a.c. power port 100 kHz repetition frequency
Surges Line-to-line	IEC 61000-4-5	±0.5 kV, ±1 kV for Input a.c. power port
Conducted disturbances induced by RF fields	IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz 6 Vrms in ISM and amateur radio bands between 150 kHz and 80 MHz 80 % AM at 1 kHz for Input a.c. power port
Rated power frequency magnetic fields	IEC 61000-4-8	30 A/m 50 Hz and 60 Hz for enclosure port

EN

11. Guidance and Manufacturer's Declaration

Phenomenon	Basic EMC standard	IMMUNITY TEST LEVELS
Voltage dips	IEC 61000-4-11	0 % U_T ; 0.5 cycle At $0^\circ, 45^\circ, 90^\circ, 135^\circ, 180^\circ, 225^\circ, 270^\circ$ and 315° for Input a.c. power port
		0 % U_T ; 1 cycle and 70 % U_T ; 25/30 cycles single phase: at 0° for Input a.c. power port
Voltage interruptions	IEC 61000-4-11	0 % U_T ; 250/300 cycle for Input a.c. power port
<p>Note: U_T is the A.C. mains voltage prior to application of the test level.</p>		

Table 3 - Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications device

Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum power (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
385	380 to 390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430 to 470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0.3	28
710	704 to 787	LTE Band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
745						
780						
810	800 to 960	GSM 800/900, TETRA 800, I den 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28
870						
930						
1720	1700 to 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
1845						
1970						
2450	2400 to 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
5240	5100 to 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9
5500						
5785						

EN

EMC tests have included the AC-adaptor as included with the product.

11. Guidance and Manufacturer's Declaration

The **Bluetooth**[®] word mark and logos are registered trademarks owned by the Bluetooth SIG, Inc. and any use of such marks by OMRON HEALTHCARE Co., Ltd. is under license. Other trademarks and trade names are those of their respective owners.

App Store and App Store logo are service marks of Apple Inc., registered in the U.S and other countries.

Google Play and the Google Play logo are trademarks of Google LLC.

12. Useful Information

What is Blood Pressure?

Blood pressure is a measure of the force of blood flowing against the walls of the arteries. Arterial blood pressure is constantly changing during the course of the heart's cycle.

The highest pressure in the cycle is called the Systolic Blood Pressure; the lowest is the Diastolic Blood Pressure. Both pressures, the Systolic and Diastolic, are necessary to enable a physician to evaluate the status of a patient's blood pressure.

EN

What is Arrhythmia?

Arrhythmia is a condition where the heartbeat rhythm is abnormal due to flaws in the bio-electrical system that drives the heartbeat. Typical symptoms are skipped heartbeats, premature contraction, an abnormally rapid (tachycardia) or slow (bradycardia) pulse.

What is Afib?

Atrial fibrillation (also called Afib or AF) is a quivering or irregular heartbeat (arrhythmia) that can lead to blood clots, stroke, heart failure and other heart-related complications. During atrial fibrillation, the heart's two upper chambers (the atria) beat chaotically and irregularly — out of coordination with the two lower chambers (the ventricles) of the heart. Episodes of atrial fibrillation can come and go, or you may develop atrial fibrillation that doesn't go away and may require treatment.

Afib indicator function detects the possibility of Afib with an accuracy of 94.2 % (with sensitivity of 95.5 % and specificity of 93.8 %) as demonstrated in the study* with Single-lead ECG as reference measurement.

* M. Ishizawa, T. Noma, T. Minamino et al. Multiple measurements with automated blood pressure monitor can detect atrial fibrillation with high sensitivity and specificity in general cardiac patients, ESC Congress 2018

Manufacturer 	OMRON HEALTHCARE Co., Ltd. 53, Kunotsubo, Terado-cho, Muko, KYOTO, 617-0002 JAPAN		
EU-representative <table border="1" data-bbox="52 768 228 827"> <tr> <td data-bbox="52 768 135 827">EC</td> <td data-bbox="135 768 228 827">REP</td> </tr> </table>	EC	REP	OMRON HEALTHCARE EUROPE B.V. Scorpius 33, 2132 LR Hoofddorp, THE NETHERLANDS www.omron-healthcare.com
EC	REP		
Importer in EU	www.omron-healthcare.com		
Production facility	OMRON HEALTHCARE MANUFACTURING VIETNAM CO., LTD. No.28 VSIP II, Street 2, Vietnam-Singapore Industrial Park II, Binh Duong Industry-Services-Urban Complex, Hoa Phu Ward, Thu Dau Mot City, Binh Duong Province, Vietnam		
Exclusive distributor and importer of OMRON medical equipment in Russia	“ComplectServise Co Ltd” 4, building 10, attic, Solnechnogorskaya street, Moscow 125413, Russia www.csmedica.ru		

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